

TECHNICAL DATA

SpectroLight 600

Laser diode	<div><input checked="" type="checkbox"/> Wavelength: 658 nm, optical power: 120 mW, adjustable</div> <div><input type="checkbox"/> other wavelengths e.g. 525 nm, 785 nm (optional)</div>
Detector	<div><input checked="" type="checkbox"/> Photomultiplier tube, dark count rate < 300 Hz quantum efficiency 5-7%, count sensitivy $1.5 \cdot 10^5$ Hz/pW</div> <div><input checked="" type="checkbox"/> For single photon counting</div> <div><input checked="" type="checkbox"/> Scattering angle 142°</div>
Correlator	<div><input checked="" type="checkbox"/> Multi-tau architecture correlator to cover a wide sample time range</div> <div><input checked="" type="checkbox"/> Sample time from 400 ns to ~ 1 s</div> <div><input checked="" type="checkbox"/> Total 208 channels, quasi logarithmic cannel spacing</div>
Sensitivity	<div><input checked="" type="checkbox"/> Sample concentration with standard laser (658 nm)</div> <div><input checked="" type="checkbox"/> Minimum 0.05 mg/ml (150 kDa), 0.1 mg/ml (30 kDa), 0.3 mg/ml (14 kDa) proteins (e.g. IgG Ab, Proteinase K, Lysozyme)</div> <div><input checked="" type="checkbox"/> Maximum concentration >150 mg/ml</div>
Imaging system	<div><input checked="" type="checkbox"/> Built-in microscope</div> <div><input checked="" type="checkbox"/> 5 magnification steps: 0.63, 1.25, 2.0, 3.2, 6.4</div> <div><input checked="" type="checkbox"/> Field of view (mm): 10.5x7.6, 5.2x2.9, 3.3x2.5, 2.0x1.5, 1.0x0.75</div> <div><input checked="" type="checkbox"/> Optical resolution per pixel: 19.5 μm, 10 μm, 6.3 μm, 3.9 μm, 2.5 μm</div> <div><input checked="" type="checkbox"/> CMOS colour camera 2048x1536 pixel</div>
Illumination	<div><input checked="" type="checkbox"/> Bright light integrated LED</div> <div><input type="checkbox"/> UV by integrated light source (optional)</div> <div><input type="checkbox"/> Colour light source (optional)</div>
Temperature control	<div><input checked="" type="checkbox"/> Built-in temperature control</div> <div><input checked="" type="checkbox"/> Temperature range 8°C to 60°C (at ambient temperature of 20°C)</div>
Sample properties	<div><input checked="" type="checkbox"/> Sample volume/well: min. ~ 20 nl, max. ~ 30 μl (depending on plate layout)</div> <div><input checked="" type="checkbox"/> Particle sizes from 0.7 nm to ~ 5 μm</div>
Sample container	<div><input checked="" type="checkbox"/> Plates in SBS format</div> <div><input checked="" type="checkbox"/> Sitting drop: e. g. MRC 96 well, Douglas Instruments Vapour Batch Plate</div> <div><input checked="" type="checkbox"/> Cryo-EM grids, microfluidic chips, customised containers</div>
Hardware	<div><input checked="" type="checkbox"/> Table top system (LxWxH) 650 mm x 270 mm x 450 mm</div> <div><input checked="" type="checkbox"/> Weight: ~ 26 kg</div> <div><input checked="" type="checkbox"/> Power consumption 115 to 230 V, 100 W</div> <div><input checked="" type="checkbox"/> Mini PC attached to the monitor (22 inch)</div>
Software features	<div><input checked="" type="checkbox"/> Operating software "SpectroCrystal" runs on Linux OpenSUSE 15.6</div> <div><input checked="" type="checkbox"/> Fully automated, drop centering, DLS positioning and plate scanning</div> <div><input checked="" type="checkbox"/> Integrated SQL database for storage and retrieval of images and DLS data</div> <div><input checked="" type="checkbox"/> Multi user workstation remote data management "Remote LIMS"</div> <div><input checked="" type="checkbox"/> Control of light source parameters</div> <div><input checked="" type="checkbox"/> Live view of the camera stream</div> <div><input checked="" type="checkbox"/> Graphical hstogramming software</div> <div><input checked="" type="checkbox"/> Radius distribution 2D and 3D</div> <div><input checked="" type="checkbox"/> Autopilot for scheduling multiple tasks</div> <div><input checked="" type="checkbox"/> Automated laser intensity admustment</div> <div><input checked="" type="checkbox"/> Diffusion interaction parameter (k_D) analysis</div> <div><input checked="" type="checkbox"/> Individual DLS scan evaluation and management</div> <div><input checked="" type="checkbox"/> Scattered Light, mass and number weithted statistics</div> <div><input checked="" type="checkbox"/> Autoscoring and user scoring options for DLS data</div> <div><input checked="" type="checkbox"/> User scoring options for image data</div> <div><input checked="" type="checkbox"/> Multiple data export funktions including customizable reporting options</div> <div><input type="checkbox"/> Integration into existing data management environment (optional)</div> <div><input type="checkbox"/> Combination with an automated plate handling system (hotel), (optional)</div>